


PED IgA ELISA

Porcine Epidemic Diarrhea virus

Porcine epidemic diarrhea virus (PEDV) is a member of the family Coronaviridae. PEDV causes acute enteritis in swine of all ages, and it is often fatal for neonatal piglets.

To protect piglets from PEDV, Sow transfers immunoglobulin through colostrums to their children until the piglets acquire adaptive immunity. And many reports suggested that IgA is important for protection of PEDV. Anigen IgA ELISA measures preventive anti PED-virus IgA titers of sow and predict defense capability of piglets induce by passive antibody transfer.



Indications


- Quantitative detection of PED IgA antibody
- Screening for defensive capacity against PEDV

Special Features

- Easy sample collection
- Optimal screening method for defensive capacity of PEDV
- Specimen : Colostrums
- Reading Time : 1 hour and 45 minutes
- Survival rate of IgA positive confirmed group after challenge : see tables below

Quick Procedure

1. Prepare PED antigen coated test plate
2. Dispense 100 μ l of sample diluents into the all strip well.
3. Dispense 10 μ l of positive, negative control and samples to each well
4. Incubate the plate at 37 \pm 1 $^{\circ}$ C for 60 minutes
5. Wash the plate at 5 times
6. Dispense 100 μ l of diluted enzyme conjugate to each well
7. Incubate the wells for 30 minutes at 37 \pm 1 $^{\circ}$ C
8. Wash the plate at 5 times
9. Dispense 100 μ l of substrate to each well
10. Incubate the wells for 15 minutes at room temperature (18~25 $^{\circ}$ C)
11. Dispense 100 μ l of stopping solution
12. Measure the optical density (OD) at 450nm with reference wavelength at 620nm
13. Cut off value = [0.35+ the mean Negative Control absorbance]



Ordering Information

Cat. No.	Description	Type	Pack size
EB44-10	PED IgA Ab ELISA	Microplate	480 Wells/Kit